

Human CD3E&CD3D (Primary Amine Labeling) protein (Recombinant) (C-His) {Biotin} (STJP020036)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Asp23-Glu120
Immunogen Sequence	DGNEEMGGITQTPYKVSISG TTVILTCPQYPGSEILWQHN DKNIGGDEDDKNIGSDEDHL SLKEFSELEQSGYYVCYPRG SKPEDANFYLYLRARVCEF K IPIEELED R V F V N C N T S I T W VEGTVGTL L S D I T R L D L G K R I L D P R G I Y R C N G T D I Y K D K E S T V Q V H Y R M C Q S C V E L D P A T V A

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---