

Human IL-23/IL-12B&IL-23A protein (Recombinant) (C-His&Avi) {Biotin} (STJP020038)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Arg20-Pro189
Immunogen Sequence	RAVPGGSSPAWTQCQQLSQK LCTLAWSAHLVGHMDLREE GDEETTNDVPHIQCGDGDGP QGLRDNSQFCLQRIHQGLIF YEKLLGSDIFTGEPSSLPDS PVGQLHASLLGLSQLLQPEG HHWETQQIPSLSPSQPWQRL LLRFKILRSLQAFVAVAAARV FIWELKKDYYVVELDWYPDA PGEMVVLTCDTPEEDGITWT LDQSSEVLGSGKTLTIQVKE FGDAGQYTCHKGGEVLSHS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---