

Human uPAR/PLAUR/CD87 protein (Recombinant) (C-His&Avi) {Biotin} (STJP020050)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	5329
Gene Symbol	PLAUR
UniProt ID	UPAR_HUMAN
Immunogen Region	Leu23-Gly305
Immunogen Sequence	LRCMQCKTNGDCRVEECALG QDLCRTTIVRLWEEGEELEL VEKSCTHSEKTNRTLSTYRTG LKITSLTEVVCGLDLCNQGN SGRAVTYSRSRYLECISCGS SDMSCERGRHQSLQCRSP EE QCLDVVTHWVQEGEEGRPKD DRHLRGCGYLPGCPGSNGFH NNDTFHFLKCCNTTKCNEGP ILELENLPQNGRQCYSCKGN STHGCSSEETFLIDCRGPMN QCLVATGTHEPKNQSYMVR
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.