

Human uPAR/PLAUR/CD87 protein (Recombinant) (C-hFc) (STJP020051)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤-70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5329
Gene Symbol	PLAUR
UniProt ID	UPAR_HUMAN
Immunogen Region	Leu23-Gly305
Immunogen Sequence	LRCMQCKTNGDCRVEECALG QDLCRTTIVRLWEEGEELEL VEKSCETHSEKTNRTLSTYRTG LKITSLTEVVCGLDLCNQGN SGRAVTYSRSRYLECI SCGS SDMSCERGRHQSLQCRSPEE QCLDVVTHWIEGEEGRPKD DRHLRGCGYLPGCPGSGNGFH NNDTFHFLKCCNTTKCNEGP ILELENLPQNGRQCYSCKGN STHGCSSEETFLIDCRGPMN QCLVATGTHEPKNQSYMVR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---