

Human C7orf50 (Cholesin) protein (Recombinant) (C-6His) (STJP020052)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84310
Gene Symbol	CHLSN
UniProt ID	CHOLN_HUMAN
Immunogen Region	Met1-Ser194
Immunogen Sequence	MAKQKRKVPVTEKKNKKLK KASAEGPLLGPEAAPSGEGA GSKGEAVLRPGLDAEPELSP EEQRVLERKLLKKERKKEERQ RLREAGLVAQHPPARRSGAE LALDYLCRWAQKHKNWRFQK TRQTWLLHMYDSKVPDEH FSTLLAYLEGLQGRARELTV QKAEALMRELDEEGSDPPLP GRAQRIRQVLQLLS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---