

Human Alkaline phosphatase (Intestinal type) /ALPI protein (Recombinant) (C-6His) (STJP020054)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	248
Gene Symbol	ALPI
UniProt ID	PPBI_HUMAN
Immunogen Region	Val20-Thr502
Immunogen Sequence	VIPAEENPAFWNRQAAEAL DAAKKLQPIQKVAKNLILFL GDGLGVPTVTATRILKGQKN GKLGPEPLAMDRFPYLALS KTYNVDRQVPDSAATATAYL CGVKANFQTIGLSAAARFNQ CNTTRGNEVISVMNRAKQAG KSVGVVTTTRVQHASPAGTY AHTVNRNWYSDADMPASARQ EGCQDIATQLISNMDIDVIL GGGRKYMFPMPGTPDPEYPAD ASQNGIRLDGKNLVQEWLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---