

Mouse CCL3/MIP-1 alpha protein (Recombinant) (C-His) (STJP020056)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO Cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in 50mM Acetate-sodium acetate pH3.6.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20302
Gene Symbol	Ccl3
UniProt ID	CCL3_MOUSE
Immunogen Region	Ala24-Ala92
Immunogen Sequence	APYGADTPTACCFSYSRKIP RQFIVDYFETSSLCSQPGVI FLTKRNRQICADSKETWVQE YITDLELNA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---