

Human Mature BMP-9/GDF-2 protein (Recombinant) (NO-Tag) (STJP020057)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2658
Gene Symbol	GDF2
UniProt ID	GDF2_HUMAN
Immunogen Region	Ser320-Arg429
Immunogen Sequence	SAGAGSHCQKTSLRVNFEDI GWDSWIIAPKEYEAYECKGG CFFPLADDVTPTKHAIVQTL VHLKFPTKVGKACCVPTKLS PISVLYKDDMGVPTLKYHYE GMSVAECGCR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---