

Human BLNK protein (Recombinant) (N-His) (STJP020061)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29760
Gene Symbol	BLNK
UniProt ID	BLNK_HUMAN
Immunogen Region	Met1-Ser456
Immunogen Sequence	MDKLNKITVPASQKLRQLQK MVHDIKNNEGGIMNKIKLK VKAPPSVPRRDYASESPADE EEQWSDDFSDYENPDEHSD SEMYVMPAEENADDSYEPPP VEQETRPVHPALPFARGEYI DNRSSQRHSPPFSKTLPSKP SWPSEKARLTSTLPALTALQ KPQVPPKPKGLLEADYVV PVEDNDENYIHPTESSPPP EKAPMVNRSTKPNSSTPASP PGTASGRNSGAWETKSPPP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---