

Human Mature BMP-10 protein (Recombinant) (NO-Tag) (STJP020063)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in 4mM HCl.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27302
Gene Symbol	BMP10
UniProt ID	BMP10_HUMAN
Immunogen Region	Ser22-Arg424
Immunogen Sequence	SPIMNLEQSPLEEDMSLFGD VFSEQDGVDFNTLLQSMKDE FLKTLNLSDIPTQDSAKVDP PEYMLELYNKFATDRTSMPS ANIIRSFKNEDLFSQPVSFN GLRKYPLLFNVSIPIHHEEVI MAELRLYTLVQRDRMIYDGV DRKITIFEVLESKGDNEGER NMLVLVSGEIYGTNSEWETF DVTDAIRRWQKSGSSTHQL E VHIESKHDEAEDASSGRLEI DTSQNKHNPLLIVFSDDQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---