

Mouse IL-17F protein (Recombinant) (C-His) (STJP020067)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	257630
Gene Symbol	Il17f
UniProt ID	IL17F_MOUSE
Immunogen Region	Met1-Ala161
Immunogen Sequence	MKCTRETAMVKSLLLMLGL AILREVAARKNPKAGVPALQ KAGNCPPLEDNTVRVDIRIF NQNQGISVPREFQNRSSSPW DYNITRDPHRFPSEIAEAQC RHSGCINAQGQEDSTMNSVA IQQEILVLRREPQGCSNSFR LEKMLLKVGCTCVKPIVHQA A

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---