

Human CCL16/HCC-4 protein (Recombinant) (C-His) (STJP020077)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4).
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6360
Gene Symbol	CCL16
UniProt ID	CCL16_HUMAN
Immunogen Region	Gln24-Gln120
Immunogen Sequence	QPKVPEWVNTPSTCCLKYYE KVLPRRLVGYRKALNCHLP AIIFVTKRNREVCTNPNDWW VQEYIKDPNLPLLPTRNLST VKIITAKNGQPQLLSQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---