

Human IL-4 protein (Recombinant) (No-tag) (STJP020080)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3565
Gene Symbol	IL4
UniProt ID	IL4_HUMAN
Immunogen Region	His25-Ser153
Immunogen Sequence	HKCDITLQEIITLNSLTEQ KTLCTELTVTDIFAASKNTT EKETFCAATVLRQFYSHHE KDTRCLGATAQQFHRHKQLI RFLKRLDRNLWGLAGLNSCP VKEANQSTLENFLERLKTIM REKYSKCSS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---