

African swine fever virus P30 protein (Recombinant) (NO-Tag) (STJP020082)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Supplied as 0.22 µg/ml filtered solution in PBS, 300 mM NaCl, 10% glycerol, (pH 7.4).
Storage Instruction	Store at -70 °C. This product is stable at ≤-70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Ile86-Lys194
Immunogen Sequence	ILHVLFEETESSASSENH EKNDNETNECTSSFELFEQ EPSSEVPKDSKLYMLAQKTV QHIEQYGGKAPDFNKVIRAHN FIQTIYGTPLKEEKEVRL MVIKLLKKK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---