

Human MMP-10 protein (Recombinant) (C-6His) (STJP020088)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of 150 mM NaCl, 5 mM CaCl ₂ , 50 mM Tris, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4319
Gene Symbol	MMP10
UniProt ID	MMP10_HUMAN
Immunogen Region	Tyr18-Cys476
Immunogen Sequence	YPLSGAAKEEDSNKDLAQY LEKYNNLEKDVKQFRRKDSN LIVKKIQGMQKFLGLEVTGK LDTDTLEVMRKPRCGVPDVG HFSSFPGMPKWRKTHLYRI VNYTPDLPRDAVDSAIEKAL KVVWEEVPLTFSRLYEGEAD IMISFAVKEHGDFYSFDGPG HSLAHAYPPGPLYGDIHFD DDEKWTEASGNTNLFVAAH ELGHSGLGFHSANTEALMYP LYNSTELAQFRLSQDDVN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---