

Human Prostatic Acid Phosphatase/ACPP protein (Recombinant) (C-6His) (STJP020091)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	0.22 Mu m filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 7.5.
Storage Instruction	Store at -70 °C. This product is stable at ≤ -70 °C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	55
Gene Symbol	ACP3
UniProt ID	PPAP_HUMAN
Immunogen Region	Lys33-Asp386
Immunogen Sequence	KELKFVTLVFRHGDRSPIDT FPTDPIKESSWPQGFGQLTQ LGMEQHYELGEYIRKRYRKFLNESYKHEQVYIRSTDVRT LMSAMTNLAALFPPEGVSIW NPILLWQPIPVHTVPLSEDQLLYLPFRNCPRFQELESETL KSEEFQKRLHPYKDFIATLG KLSGLHGQDLFGIWSKVYDPLYCESVHNFTLPSWATEDTM TKLRELSELSLLSLYGIHKQ KEKSRLQGGVLVNEILNHM
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.