

Human CMRF35-like molecule 5/CD300LD protein (Recombinant) (C-hFC) (STJP020094)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100131439
Gene Symbol	CD300LD
UniProt ID	CLM5_HUMAN
Immunogen Region	Ala19-His165
Immunogen Sequence	AKITGPTTVNGSEQGSLTVQ CAYGSGWETYLKWRCQGADW NYCNILVKTNNGSEQEVKKNR VSIRDNQKNHVFTVTMENLK RDDADSYWCGTERPGIDLV KVQVTINPGTQTAVSEWTTT TASLAFTAAATQKTSSPLTR SPLKSTH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---