

Human STAT1 protein (Recombinant) (NO-Tag) (STJP020096)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, 2mM DTT 1mM EDTA pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6772
Gene Symbol	STAT1
UniProt ID	STAT1_HUMAN
Immunogen Region	Met1-Val712
Immunogen Sequence	MSQWYELQQLDSKFLEQVHQ LYDDSPMEIRQYLAQWLEK QDWEHAANDVSFATIRFHD LSQLDDQYSRFSLENNFLLQ HNIRKSKRNLDNFQEDPIQ MSIIYSCLKEERKILENAQ RFNQAQSGNIQSTVMLDKQK ELDSKVRNVKDKVMCIEHEI KSLEDLQDEYDFKCKTLQNR EHETNGVAKSDQKQEQLLLK KMYLMLDNKRKEVVKIIEI LNVTELTQNALINDELVEW

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---