

Human N-Cadherin/CDH2 protein (Recombinant) (C-hFC) (STJP020098)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1000
Gene Symbol	CDH2
UniProt ID	CADH2_HUMAN
Immunogen Region	Asp160-Ala724
Immunogen Sequence	DWVIPPINLPENSRGPPFQE LVRI RSDRDKNLSLRYSVTG PGADQPPTGIFIINPISGQL SVTKPLDREQIARFHLRAHA VDINGNQVENPIDIVINVID MNDNRPEFLHQVWNGTVPEG SKPGTYVMTVTAIDADDPNA LNGMLRYRIVSQAPSTPSPN MFTINNETGDIITVAAGLDR EKVQYQYTLIIQATDMEGNPT YGLSNTATAVITVDVNDNP PEFTAMTFYGEVPENRVDI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---