

## Human FCRL2/FCRH2/IRTA4 protein (Recombinant) (C-6His) (STJP020099)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">79368</a>
Gene Symbol	<a href="#">FCRL2</a>
UniProt ID	<a href="#">FCRL2_HUMAN</a>
Immunogen Region	Leu20-Asp395
Immunogen Sequence	LTLVAPSSVFEGDSIVLKCQ GEQNWKIQKMAYHKDNKELS VFKKFSDFLIQSAVLSDSGN YFCSTKGQLFLWDKTSNIVK IKVQELFQRPVLTASSFQPI EGGPVSLKCETRLSPQLDV QLQFCFFRENQVLGSGWSSS PELQISAVWSEDTGSYWCKA ETVTHRIRKQSLQSQIHVQR IPISNVLSLEIRAPGGQVTEG QKLILLCSVAGGTGNVTFWSW YREATGTSMGKKTQRSLSA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---