

Human CCL14/HCC-1/HCC-3 (28-93) protein (Recombinant) (C-6His) (STJP020100)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6358
Gene Symbol	CCL14
UniProt ID	CCL14_HUMAN
Immunogen Region	Gly28-Asn93
Immunogen Sequence	GPYHPSECCFTYTTYKIPRQ RIMDYETNSQCSKPGIVFI TKGHSVCTNP SDKWVQDYI KDMKEN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---