

Human CR2/CD21 protein (Recombinant) (C-Avi&6His) (STJP020102)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1380
Gene Symbol	CR2
UniProt ID	CR2_HUMAN
Immunogen Region	Ile21-Arg971
Immunogen Sequence	ISCGSPPPILNGRISYYSTP IAVGTVIRYSCSGTFRIDGE KSLLCITKDKVDGTWDKPAP KCEYFNKYSSCPEIIVPGGY KIRGSTPYRHGDSVTFACKT NFSMNGNKS VWCQANNMWGP TRLPTCVSVFPLECPALPMI HNGHHTSENVGSIAPGLSVT YSCESGYLLVGEKIINCLSS GKWSAVPPTCEEARCKSLGR FPNGKVKEPPILRVGVTANF FCDEGYRLQGPPSSRCVIA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---