

Human Alkaline phosphatase/tissue-nonspecific isozyme/ALPL protein (Recombinant) (C-6His) (STJP020103)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.2 Mu m filtered solution of 20mM HEPES, 150mM NaCl, 2mM MgSO4, 0.1mM ZnCl2, pH 7.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	249
Gene Symbol	ALPL
UniProt ID	PPBT_HUMAN
Immunogen Region	Leu18-Ala500
Immunogen Sequence	LVPEKEKDPKYWRDQAQETL KYALELQKLNTNVAKNVIMF LGDGMGVSTVTAAIRILKGQL HHNPGEETRLEMDKFPFVAL SKTYNTNAQVPDSAGTATAY LCGVKANEGTVGVSAATERS RCNTTQGNEVTSILRWAKDA GKSVGIVTTTRVNHATPSAA YAHSADRDWYSDNEMPPEAL SQGCKDIAYQLMHNIRDIDV IMGGGRKMYMPKNKTDVEYE SDEKARGTRLDGLDLVDTW
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.