

Human Insulin Receptor/INSR/CD220 protein (Recombinant) (C-Avi&His) (STJP020104)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	CHO Cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	3643
Gene Symbol	INSR
UniProt ID	INSR_HUMAN
Immunogen Region	His28-Val944
Immunogen Sequence	HLYPGEVCPGMDIRNNLTRL HELENC SVIEGHLQILLMFK TRPEDFRDLSFPKLIMITDY LLLFRVYGLES LKDLFPNLT VIRGSRLFFNYALVIFEMVH LKELGLYNLMNITRGSVRIE KNNELCYLATIDWSRILDSV EDNYIVLNKDDNEECGDICP GTAKGKTNCPATVINGQFVE RCWTHSHCQKVCPTICKSHG CTAEGLCCHSECLGNCSQPD DPTKCVACRNFYLDGRCVE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---