

Human DLL3 protein (Recombinant) (N-His) (STJP020109)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 200 mM Arginine, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10683
Gene Symbol	DLL3
UniProt ID	DLL3_HUMAN
Immunogen Region	Ala27-Arg491
Immunogen Sequence	AGVFELQIHSFGPGPGPGAP RSPCSARLPCRLFFRVCLKP GLSEEAESPICALGAALSAR GPVYTEQPGAPAPDLPLPDG LLQVPFRDAWPGTFFSFIET WREELGDQIGGPAWSLLARV AGRRRLAAGGPWARDIQRAG AWELRFSYRARCPEPPAVGTA CTRLCRPRSAPSRGCGPLRP CAPLEDECEAPLVCRAGCSP EHGFCCEQPGECRCLEGWTGP LCTVPVSTSSCLSPRGPS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------