

Human B7-H1/PD-L1/CD274 protein (Recombinant) (C-mFc) (STJP020111)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 29126 |
| Gene Symbol | CD274 |
| UniProt ID | PD1L1_HUMAN |
| Immunogen Region | Phe19-Arg238 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|