

Human Nectin-4/PVRL4 protein (Recombinant) (C-hFc) (STJP020112)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81607
Gene Symbol	NECTIN4
UniProt ID	NECT4_HUMAN
Immunogen Region	Gly32-Ser349
Immunogen Sequence	GELETSDVVTVVLGQDAKLP CFYRGDSGEGVGVAVARVD AGEQAQELALLHSKYGLHVS PAYEGRVEQPPPPRNPLDGS VLLRNAVQADEGEYECRVST FPAGSFQARLRRLRVLPPLP SLNPGPALEEGQLTLAASC TAEGSPAPSVTWDEVKGTT SSRSFKHSRSAAVTSEFHLV PSRSMNGQPLTCVVSHPGLL QDQRITHILHVSFLAEASVR GLEDQNLWHIGREGAMLKC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---