

Cynomolgus CD3E&CD3D protein (Recombinant) (C-hFc) (STJP020114)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Asp20-Arg687
Immunogen Sequence	QDGNEEMGSITQTPYQVSIS GTTVILTCSQHLGSEAQWQH NGKNKEDSGDRLFLPEFSEM EQSGYYVCYPRGSNPEDASH HLYLKARVCENCMEMDFKIP VEELED RVFKCNTSVTWVE GTVGTLTNNTRLDLGKRIL DPRGIYRCNGTDIYKD KESA VQVHYRMCQNCVELDPATLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---