

SARS-CoV-2 Spike S1 (B.1.1.529/Omicron) protein (Recombinant) (C-His&Avi) {Biotin} (STJP020128)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Gln14-Arg683
Immunogen Sequence	Gln14-Arg683 (A67V, H V69-70del, T95I, G14 2D, VYY143-145del, N 211del, L212I, ins21 4EPE, G339D, S371L, S373P, S375F, K417N, N440K, G446S, S477N, T478K, E484A, Q493 R, G496S, Q498R, N50 1Y, Y505H, T547K, D6 14G, H655Y, N679K, P 681H)

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---