

## SARS-CoV-2 Spike S1+S2 ECD (S-ECD) (B.1.1.529/Omicron) Trimer protein (Recombinant) (C-His) (STJP020129)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	Val16-Glu1213
Immunogen Sequence	Val16-Glu1213 (A67V, HV69-70del, T95I, G1 42D, VYY143-145del, N211del, L212I, ins2 14EPE, G339D, S371L, S373P, S375F, K417N, N440K, G446S, S477 N, T478K, E484A, Q49 3R, G496S, Q498R, N5 01Y, Y505H, T547K, D 614G, H655Y, N679K, P681H, N764K, D796Y, N8

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---