

Human CXCL16/SR-PSOX protein (Recombinant) (C-His) (STJP020131)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58191
Gene Symbol	CXCL16
UniProt ID	CXL16_HUMAN
Immunogen Region	Met1-Thr205
Immunogen Sequence	MGRDLRPGSRVLLLLLLLL VYLTQPGNGNEGSVTGSCYC GKRISSDSPPSVQFMNRLRK HLRAYHRCLYYTRFQLLSWS VCGGNKDPWVQELMSCLDLK ECGHAYSGIVAHQKHLLPTS PPISQASEGASSDIHTPAQM LLSTLQSTQRPTLPVGSLS DKELTRPNETTIHTAGHSLA AGPEAGENQKQPEKNAGPTA RTSAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---