

Human GFR alpha-like/GFRAL protein (Recombinant) (C-His&Avi) (STJP020134)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22um filtered solution in 20mM PB, 300mM NaCl (pH 6.5). Normally 5 % trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	389400
Gene Symbol	GFRAL
UniProt ID	GFRAL_HUMAN
Immunogen Region	Ser19-Glu351

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---