

Mouse CD3E&CD3D protein (Recombinant) (C-hFc) (STJP020138)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from 0.22um filtered solution in PBS (pH 7.4). Normally 5 % trehalose is added as protectant before lyophilization.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12500 12501
Gene Symbol	Cd3d Cd3e
UniProt ID	CD3D_MOUSE CD3E_MOUSE
Immunogen Region	Asp23-Asp108
Immunogen Sequence	DAENIEYKVSISGTSVELTC PLSDENLKWEKNGOELPQK HDKHLVLQDFSEVEDSGYYV CYTPASNKNTYLYLKARVCE YCVEVDFKIQVTEYEDKVFV TCNTSVMHLDGTVEGWFAKN KTLNLGKGVLDPRGIYLCNG TEQLAKVVSSVQVHYRMCQN CVELDSGTMA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---