

Human IL-17RA/CD217 protein (Recombinant) (C-hFc&Avi) {Biotin} (STJP020176)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	23765
Gene Symbol	IL17RA
UniProt ID	I17RA_HUMAN
Immunogen Region	Leu33-Trp320
Immunogen Sequence	LRLDHRALVCSQPGLNCTV KNSTCLDDSWIHPRNLTSS PKDLQIQLHFAHTQQGDLFP VAHIEWTLQTDASILYLEGA ELSVLQLNTNERLCVRFEFL SKLRHHHRRWRFTFSHFVVD PDQEYEVTVHHLPKIPDGD PNHQSKNFLVPDCEHARMKV TTPCMSSGSLWDPNITVETL EAHQLRVSF TLWNESTHYQI LLTSFPHMENHSCFEHMHHI PAPRPEEFHQRSNVTLLR
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.