

Human Integrin alpha V beta 5 (ITGAV&ITGB5) protein (Recombinant) (C-His&Avi) (STJP020186)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3685 3693
Gene Symbol	ITGAV ITGB5
UniProt ID	ITAV_HUMAN ITB5_HUMAN
Immunogen Region	Phe31-Val992

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---