

Cynomolgus Mesothelin/MSLN protein (Recombinant) (C-His) (STJP020191)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Asp296-Gly580
Immunogen Sequence	DVERTTCPPPEKEVHEIDENL IFYKKRELEACVDAALLAAQ MDRVDAIPFTYEQLDVLKHK LDELYPQGYPESVIRHLGHL FLKMSPEDIRKWNVTSLETL KALLKVSKGHEMSAQVATLI DHVVVGRGQLDKDTADTLTA FCPGCLCSLSPERLSSVPPS VIGAVRPQDLDTGPRQLDV LYPKARLAFQNMMSGSEYFVK IRPFLGGAPTEDLKALSQQN VSMDLATFMKLRREAVLPL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---