

Human Mesothelin/MSLN protein (Recombinant) (C-hFc&Avi) {Biotin} (STJP020194)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10232
Gene Symbol	MSLN
UniProt ID	MSLN_HUMAN
Immunogen Region	Glu296-Gly580
Immunogen Sequence	EVEKTACPSGKKAREIDESL IFYKKWELEACVDAALLATQ MDRVNAIPFTYEQLDVLKHK LDELYPQGYPESVIQHLGYL FLKMSPEDIRKWNVTSLLETL KALLEVNKGHEMSPQVATLI DRFVKGRGQLDKDTLDTLTA FYPGYLCSLSPHEELSSVPPS SIWAVRPQDLDTCDPRQLDV LYPKARLAFQNMNGSEYFVK IQSFLGGAPTEDLKALSQQN VSMDLATFMKLRRTDAVLPL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---