

Human NKG2A&KLRD1/CD159a&CD94 protein (Recombinant) (N-His&Avi) (STJP020196)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3821 3824
Gene Symbol	KLRC1 KLRD1
UniProt ID	NKG2A_HUMAN KLRD1_HUMAN
Immunogen Region	Ser34-Ile179

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---