

Human FGFR-2 alpha (IIIb) /KGFR/CD332 protein (Recombinant) (C-hFc) (STJP020217)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
UniProt ID	FGFR2_HUMAN
Immunogen Region	Arg22-Glu378
Immunogen Sequence	RPSFSLVEDTTLEPEEPPTK YQISQPEVYVAAPGESLEVR CLLKDAAVISWTKDGVHLGP NNRTVLIGEYLQIKGATPRD SGLYACTASRTVDSETWYFM VNVTDAISSGDDDDTDGAE DFVSENSNNKRAPYWTNTEK MEKRLHAVPAANTVKFRCPA GGNPMPTMRWLKNGKEFKQE HRIGGYKVRNQHWSLIMESV VPSDKGNYTCVVENEYGSIN HTYHLDVVERSPPHRPILQA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---