

Human FGFR-2 beta (IIlb) /KGFR/CD332 protein (Recombinant) (C-His) (STJP020218)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
UniProt ID	FGFR2_HUMAN
Immunogen Region	Arg152-Glu378
Immunogen Sequence	RAPYWTNTEKMEKRLHAVPA ANTVKFRCPAGGNPMPTMRW LKNGKEFKQEHRIGGYKVRN QHWSLIMESVVP SDKGNYTC VVENEYGSINHTYHLDVVER SPHRPILQAGLPANASTVVG GDVEFVCKVYSDAQPHIQWI KHVEKNGSKYGPDGLPYLKV LKHSGINSSNAEVLALFNVT EADAGEYICKVSNYIGQANQ SAWLTVLPKQQAPGREKEIT ASPDYLE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---