

Human Siglec-10 protein (Recombinant) (C-His&Avi) {Biotin} (STJP020239)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20 mM Tirs 150 mM NaCl, 200 mM L-arginine, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89790
Gene Symbol	SIGLEC10
UniProt ID	SIG10_HUMAN
Immunogen Region	R119A

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---