

## Human TGF-beta receptor type-1/TGFR-1/ALK-5 protein (Recombinant) (C-mFc&Avi) {Biotin} (STJP020244)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7046</a>
Gene Symbol	<a href="#">TGFBRR1</a>
UniProt ID	<a href="#">TGFR1_HUMAN</a>
Immunogen Region	Leu34-Glu125

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---