

Human B7-H7/HHLA2 protein (Recombinant) (C-hFc) (STJP020257)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	11148
Gene Symbol	HHLA2
UniProt ID	HHLA2_HUMAN
Immunogen Region	Ile23-Asn344
Immunogen Sequence	IFPLAFFIYVPMNEQIVIGR LDEDIILPSSFERGSEVVIH WKYQDSYKVVHSYYKGS DHLE SQDPRYANRTSLFYNEIQNG NASLFFRRVSLLDDEGIYTCY VGTAIQVITNKVVLKVG VFL TPVMKYEKRNTNSFLICSVL SVYPRPIITWKMDNTPISEN NMEETGSLDSFSINSPLNIT GSNSSYECTIENSLLKQWTW GRWTMKDGLHKMQSEHVSLSCQPVNDYFSPNQDFKVTWS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---