

Human IL-17Rc protein (Recombinant) (C-hFc) (STJP020259)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	84818
Gene Symbol	IL17RC
UniProt ID	I17RC_HUMAN
Immunogen Region	Leu21-Arg467
Immunogen Sequence	LERLVGPQDATHCSPGLSCR LWSDILCLPGDIVPAGPV LAPTHLQTELVLRCQKETDC DLCLRVAVHLAVHGHWEEPE DEEFKFGGAADSGVEEPRNAS LQAQVLSFQAYPTARCVLL EVQVPAALVQFGQSVGSVYVY DCFEAALGSEVRIWSYQPR YEKELNHTQQLPDCRGLEVW NSIPSCWALPWLNV SADGN VHLVLNVSEEQHFGLSLYWN QVQGPPKPRWHKNLTGPQI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---