

Human GFR alpha-like/GFRAL protein (Recombinant) (C-hFc) (STJP020260)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	389400
Gene Symbol	GFRAL
UniProt ID	GFRAL_HUMAN
Immunogen Region	Ser19-Glu351
Immunogen Sequence	SQTNNCTYLREQCLRDANGC KHAWRVMEDACNDSDPGDPC KMRNSSYCNLSIQYLVESNF QFKECLCTDDFYCTVNKLLG KKCINKSDNVKEDKFKWNLT TRSHHGFKGMWSCLEVAEAC VGDVVCNAQLASYLKACSAN GNPCDLKQCQAAIRFFYQNI PFNIAQMLAFCDCAQSDIPC QQSKEALHSKTCVNMVPPP TCLSVIRSCQNDELCCRHYR TFQSKCWQRVTRKCHEDEN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---