

Cynomolgus Fc gamma RIII/CD16 protein (Recombinant) (C-His) (STJP020274)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Immunogen Region	Gly17-Gln208
Immunogen Sequence	GMRAEDLPKAVVFLEPQWYR VLEKDRVTLKCGAYSPEDN STRWFHNESLISSQTSSYFI AAARVNNSGEYRCQTSLSLTL SDPVQLEVHIGWLLQAPRW VFKEEESIHLRCHSWKNTLL HKVTYLQNGKGRKYFHQNSD FYIPKATLKDSGSYFCRGLI GSKNVSSETVNITITQDLAV SSISFFPPGYQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---