

Cynomolgus NKp46/NCR1/CD335 protein (Recombinant) (C-His) (STJP020276)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
|---------------------|--|

TARGET INFORMATION

| | |
|------------------|--------------|
| Immunogen Region | Ala22-Asn254 |
|------------------|--------------|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|