

Cynomolgus/Rhesus macaque TNFRSF9/4-1BB protein (Recombinant) (C-His) (STJP020277)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Immunogen Region	Leu24-Gln186
Immunogen Sequence	LQDLCSNCPAGTFCDNNRSQ ICSPCPPNSFSSAGGQRTCD ICRQCKGVFKTRKECSSTSN AECDCISGYHCLGAECMCE QDCKQGQELTKKGCKDCCFG TFNDQKRGICRPWTNCSLDG KSVLVNGTKERDVVCGPSFA DLSPGASSATPPAPAREPGH SPQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---