

Cynomolgus Glypican-3/GPC3 protein (Recombinant) (C-His) (STJP020282)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Immunogen Region	Gln25-His559
Immunogen Sequence	QPPPPPPDATCHQVRSFFQR LQPGLKWWPETPVPGSDLQV CLPKGPTCCSRKMEEKYQLT ARLNMEQLLQSASMEKFLI IQNAAVFQEA FEIVVRHAKN YTNAMFKNNYPSLTPQAFEF VGEFFTDVSLYILGSDINVD DMVNELFDSLFPVIYTQLMN PGLPDSALDINECLRGARRD LKVFGNFPKLIMTQVSKSLQ VTRIFLQALNLGIEVINTTD HLKFSKDCGRMLTRMWWYCS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---