

Human B7-H7/HLA2 protein (Recombinant) (C-His) (STJP020301)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
---------------------	--

TARGET INFORMATION

Gene ID	11148
Gene Symbol	HLA2
UniProt ID	HLA2_HUMAN
Immunogen Region	Ile23-Asn344

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---